




IN THE EVENT THAT FIXTURE LOCATIONS NEED TO CHANGE, PLEASE CHECK WITH SKYTRON ON THE REVISED LAYOUT LOCATION. FINAL CONSTRUCTION DOCUMENTS SHOULD BE PROVIDED BY THE ARCHITECT.



**\*GENERIC LAYOUT\***

EQUIPMENT SCHEDULE	
ID	DESCRIPTION
A	EQUIPMENT / DUAL DISPLAY / LIGHT
* IF SKYVISION IS APPLICABLE REFER TO THE LAST PORTION OF BOOM PACKAGE*	




INITIAL: \_\_\_\_\_  
DATE: \_\_\_\_\_

PROJECT #: 21-000  
PRELIMINARY  
PLOT DATE: 9/08/2021

ENDOSCOPY / GI

ROOM TYPE: ENDO  
DESCRIPTION: ROOM LAYOUT  
REV #: 0

SHEET  
**00**



REFER TO AUA3 DATA SHEET  
FOR SURGICAL LIGHT  
SPECIFICATIONS



INITIAL: \_\_\_\_\_  
DATE: \_\_\_\_\_

**Max Fixture Weight:** 1,188 lbs.  
**Max Moment Load:** 6,922 ft. lbs.

**Equipment Capacity:** 4CMC-49 - 351 lbs.  
4FSB-2 - 83 lbs.

**ENDOSCOPY / GI**



REF #: C-75160  
MODEL: F320 SERIES  
DESCRIPTION: ELEVATION DETAILS

QTY: 1  
REV #: 0  
SHEET  
**A1**

FRONT-LEFT VIEW

FRONT-RIGHT VIEW

Note:  
Communications can only  
be placed above the  
top-most separator plate



DESCRIPTION	QTY.
OXYGEN	1
MEDICAL AIR	1
VACUUM	2
CARBON DIOXIDE (DISS)	1

ELECTRICAL OUTLETS: (8) 125V, 20A DUPLEX - RED

GAS FACEPLATE STYLE: DISS

CARRIER DIMENSIONS: 56"H x 22"W x 28.75"D

GAS COLOR KEY:

CARBON DIOXIDE	MEDICAL AIR	OXYGEN
HELIOX	NITROGEN	VACUUM
INSTRUMENT AIR	NITROUS OXIDE	WAGD

INITIAL: \_\_\_\_\_  
DATE: \_\_\_\_\_

FOR A COMPLETE LIST OF ACCESSORIES  
REFER TO YOUR QUOTE

ENDOSCOPY / GI

QTY: 1  
REV #: 0  
REF #: C-75160  
MODEL: F320 SERIES  
DESCRIPTION: CARRIER DETAILS

SHEET  
A2